

# INStest™

## SARS-COV-2 RT-qPCR Assay For Confident COVID-19 Diagnosis

### Triple Fluorescence RT-PCR Method

- INSTEST SARS-COV-2 RT-qPCR Assay is a real-time (rt) reverse transcriptase (RT) polymerase chain reaction (PCR) test intended for the qualitative detection of ORF1ab and N genes of the SARS-CoV-2 in nasopharyngeal (NP) swab, oropharyngeal (OP) swab and sputum specimen. These specimens should be collected by a healthcare professional after a medical provider advises the patient to undergo testing.
- The SARS-COV-2 RT-qPCR Assay is applicable to fluorescent quantitative PCR Analyzers with FAM, HEX/VIC, TEXAS, and RED/ROX channels.

#### I High Accuracy

Limit of Detection: 200 copies/mL  
Sensitivity: 100 %  
Specificity: 99.1%  
Accuracy: 99.3%

#### I Convenient Operation

Triplex-detections in a single tube  
(ORF1ab & N genes, human RNase P (RNP) gene)  
Primer and probe premixed liquid form

#### I Reliable

Anti-pollution system: UDG  
Internal control: Human RNaseP (RNP) gene  
Positive control: Pseudovirus

#### I Fast

Result in about 1 hour





### 1. Specimen Collection

**Product:** Virus Specimen Stabilizer

**Cat.No.:** VS-ACOP50

**Storage and Stability:**

15 ~ 25°C for 12 months.

#### Advantages:

Quickly inactivate the virus, eliminate the risk of infection

Protects the integrity of the nucleic acid



### 3. PCR Testing

**Product:** SARS-COV-2 RT-qPCR Assay

**Cat.No.:** ICO-ACOP50

**Limit of Detection:** 200 copies/mL

**Storage and Stability:** stored at -30 ~ -15°C with protection from light for 6 months.

\* The test can eliminate the non-specific interference of SARS2003 strain and BatSARS-like virus strain, ensuring accurate detection of SARS-COV-2 virus.



### 2. Sample Processing (Virus RNA Extraction)

**Product:**

Virus Nucleic Acid Extraction Kit (spin column)

**Cat.No.:** VNE-ACOP50

**Storage and Stability:**


Room Temperature for 12 months.

#### Advantages:


**Fast:** Obtain high-purity nucleic acid within 10 minutes

| Components                | Amount 50 tests | Amount 100 tests |
|---------------------------|-----------------|------------------|
| Lysis Solution            | 18 mL           | 36 mL            |
| Washing Buffer            | 15 mL           | 30 mL            |
| Elution Buffer            | 3 mL            | 6 mL             |
| Adsorption Columns        | 50 tubes        | 100 tubes        |
| Collection Tubes (2 mL)   | 50 tubes        | 100 tubes        |
| Collection Tubes (1.5 mL) | 50 tubes        | 2 x 50 tubes     |
| Package Insert            | 1               | 1                |

| Components       | Amount 50 tests | Amount 100 tests | Ingredient  | Cap Color |
|------------------|-----------------|------------------|---|-----------|
| Reaction Mix     | 675μl*2 tubes   | 900μl*3 tubes    | Orf1ab, N genes and Human RNaseP (RNP) gene primers and probes, Buffer, deoxy-ribonucleoside triphosphate (dNTPs) | Red       |
| Enzyme Mix       | 200μl*1 tube    | 400μl*1 tube     | RNase Inhibitor, UDG, Reverse Transcriptase, Taq DNA polymerase   | Blue      |
| Positive Control | 125μl*1tube     | 250μl*1tube      | RAN Pseudovirus Containing Target Gene (Orf1ab, N genes and Human RNaseP (RNP) gene )                             | Yellow    |
| Negative Control | 125μl*1tube     | 250μl*1tube      | DEPC-Treated Water  | Green     |
| Package insert   | 1               | 1                | /   | /         |



Test results of positive sample



Test results of weak positive sample

### Ordering Information

| Cat. No.   | Product Description               | Pack      | CE Status |
|------------|-----------------------------------|-----------|-----------|
| ICO-ACOP50 | SARS-COV-2 RT-qPCR Assay          | 50 T/100T | CE        |
| VNE-ACOP50 | Virus Nucleic Acid Extraction Kit | 50 T/100T | CE        |
| VS-ACOP50  | Virus Specimen Stabilizer         | 50 T/100T | CE        |